

## FHL Transfer

### SOFTFIX-PK<sup>®</sup> Interference Screw



# Technique Guide



## Step 1

Make a 5-7cm longitudinal incision just medial to the Achilles tendon. Deeper dissection is carried out through the fascia until the FHL tendon and its muscle belly is identified. Movements of the great toe identify the FHL tendon.

The FHL tendon is dissected out from the surrounding soft tissues as distally as possible and released.

## Step 2

The tendon is stitched using #2 BioFiber<sup>®</sup> or BioFiber<sup>®</sup> Loop.



## Step 3


Check the diameter of the tendon using a graft sizer for foot and ankle.

## Step 4

Place the 2.4mm Drill Tip Suture Passing Pin into the dorsal medial aspect of the calcaneus.

**Note:** If the tendon diameter is less than 4.0mm, use a 1.2mm threaded guidewire to drill calcaneus.





## Step 5


Over ream, the Guide Pin with a fluted cannulated reamer until the desired depth.

**Note:** On creation of bone tunnel,

- (1) Proper screw and drill hole diameter depend on the exact diameter of the tendon graft.
- (2) Bone tunnel diameter should be 1.0mm larger than the tendon's size (e.g. 5.0mm graft requires 6.0mm diameter tunnel).
- (3) Bone tunnel depth should be 2.0mm longer than the length of the screw (6mm x 23mm screw = depth of 25mm).

## Step 6

Plantar flex the foot and mark the tendon at the insertion site while under tension. Put a mark over the FHL tendon depending on the tunnel depth. Pass FHL graft in the tunnel and pull it up to the mark.



## Step 7

Insert 1.2mm guidewire, Blunt Tip (Cat. No.10096) for dia. 5.0mm-8.0mm SOFTFIX-PK® interference screws into the bone tunnel.

**Note:** Use 0.8mm guidewire, Blunt Tip (Cat. No.10094) for dia. 4.0mm SOFTFIX-PK® interference screw.



## Step 8

Insert the appropriately sized SOFTFIX-PK® Interference screw over the guidewire and advance it into the bone tunnel using an appropriate screwdriver. The screw is seated precisely only when it is flushed inside the cortical bone. Remove the screwdriver and guidewire. Cut the remaining suture.

## Ordering Information

| Catalog No. | Product Description  |
|-------------|--|
| 10053       | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 4.0mm, Len. 10.0mm |
| 10054       | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 4.0mm, Len. 15.0mm |
| BAK-7164.15 | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 5.0mm, Len. 15mm   |
| BAK-7164.20 | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 5.0mm, Len. 20mm   |
| BAK-7165.23 | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 6.0mm, Len. 23mm   |
| BAK-7166.23 | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 7.0mm, Len. 23mm   |
| BAK-7167.23 | SOFTFIX-PK® Interference screw, cannulated, roundhead; Dia. 8.0mm, Len. 23mm   |

## Related Instrumentation

| Catalog No. | Product Description   |
|-------------|---|
| 10094       | Guide Wire, Blunt Tip for Dia. 4.0mm SOFTFIX-PK® Interference Screws Dia. 0.8mm; Len. 200mm (SS)            |
| 10096       | Guide Wire, Blunt Tip for Dia. 5.0mm-8.0mm SOFTFIX-PK® Interference Screws Dia. 1.2mm; Len. 200mm (NITINOL) |
| 10007       | Screwdriver for SOFTFIX-PK® Interference Screw Len. 23.0mm, Cannulated                                      |
| 10008       | Screwdriver for SOFTFIX-PK® Interference Screw, Cannulated, Dia. 4.0mm                                      |
| BAK-7223    | Screwdriver for SOFTFIX-PK®, Dia: 5.0mm Interference Screws, Cannulated                                     |
| BAL-024     | Threaded Guidewire Dia. 1.2; Len. 180mm, SS   |
| 10057       | Graft Sizer for Foot & Ankle, Dia. 4-5-6-7-8 mm   |
| 10085       | Fluted Reamer Cannulated, Dia. 4.0mm  |
| 10086       | Fluted Reamer Cannulated, Dia. 5.0mm  |
| 10087       | Fluted Reamer Cannulated, Dia. 6.0mm  |
| 10088       | Fluted Reamer Cannulated, Dia. 7.0mm  |
| 10089       | Fluted Reamer Cannulated, Dia. 8.0mm  |



**BIOTEK**<sup>®</sup>  
we put life in metal



**BIOTEK - Chetan Meditech Pvt. Ltd.**

Plot No: MD 04, Near Teva Company, Sanand GIDC Gate No. 2,  
Charal Gaam, Ahmedabad-382110. Gujarat, INDIA.  
Email: [contact@biotekortho.com](mailto:contact@biotekortho.com) Website: [www.biotekortho.com](http://www.biotekortho.com)

An ISO 13485 : 2016 Company  
All Implants specified in the catalogue are CE certified

