en POLYLACTIC INTERFERENCE SCREW
Instructions for use

DEVICE DESCRIPTION

The polyactic interference screw is a fixation system designed to ensure fixation of a transplant within a bone tunnel as part of a knee ligament reconstruction procedure. It is a resorbable device composed of poly(L-Lactide-co-D-Lactide) acid. The polyactic interference screw is available in different diameters (D) and lengths (L) allowing the adjustment of the size of the screw to the size of the bone tunnel and therefore an optimal fixation of the transplant. Its internal print is unique and specific to the installation instrumentation to be used for its implantation.

COMPOSITION

Component	Proportion
Poly (L-Lactide-co-D-Lactide) acid	100%

INTENDED USE

The polyactic interference screw is intended to be used in knee ligament reconstruction procedures. It ensures the fixation of the transplant within the bone tunnel during the healing period and for a maximum duration of 30 months.

SURGICAL INDICATIONS

The polyactic interference screw is indicated in the treatment of Anterior Cruciate Ligament (ACL) reconstruction.

TARGET POPULATION

The polyactic interference screw is exclusively reserved for adults, except pregnant or breastfeeding women.

The use of the polyactic interference screw has not been validated for use in the pediatric population or on pregnant or breastfeeding women.

CONTRAINDICATIONS

The use of the device is contraindicated in the following cases:

- Acute infections
- Known allergies to the components (if an allergy of any kind is suspected, relevant preliminary medical examinations must be carried out)
- Poor or inadequate bone quality (in particular in case of tumors and severe osteoporosis)
- Limited blood flow at the grafting site (necrosis, previous infections)
- Growing patients with cartilage activity
- Fever
- Local inflammation

ADVERSE EVENTS

Device related adverse events:

- Inflammation
- Cysts
- Allergic reaction
- Breakage
- Infection
- Meniscal tears
- Graft failure/instability

Procedure related adverse events:

- Arthrofibrosis/Decreased patient mobility
- Synovial damage
- Hematoma
- Bursitis
- Patellar pain syndrome
- Traumatic shock
- Tendon disorder
- Muscle strain
- Ligament rupture

Some of adverse events may be treated by drugs or require surgical reintervention.

EXPECTED CLINICAL BENEFITS

The polyactic interference screw supports ACL knee ligament reconstruction by acting as a fixation system to hold the transplant in the bone tunnel, stabilizing the knee and allowing the patient to resume activities.

STERILISATION

The polyactic interference screw is sterilized by beta irradiation in its packaging. To avoid any modification of its performances, it should never be subjected to a second sterilization before use.

CONSERVATION

The polyactic interference screw must be stored in its original undamaged packaging and kept away from moisture at a temperature lower than 32°C (90°F).

WARNINGS

DO NOT USE THE PRODUCT IF:

- PACKAGING IS OPENED OR DAMAGED
- EXPIRY DATE IS EXCEEDED.

The polyactic interference screw is a single use implantable medical device. It is for a single patient and must never be re-used. The re-use of the device may alter its performances and may lead to postoperative infections. The polyactic interference screw is a single use implantable medical device. It is for a single patient and must never be re-used. The re-use of the device may alter its performances and may lead to postoperative infections. The polyactic interference screw is a single use implantable medical device. It is for a single patient and must never be re-used. The re-use of the device may alter its performances and may lead to postoperative infections. The polyactic interference screw must only be used in an operating room.

The polyactic interference screw cannot be observed by X-ray. It does not create magnetic resonance imaging artifacts and therefore can be used safely in a magnetic resonance environment.

PRECAUTIONS FOR USE

The polyactic interference screw must be implanted by a surgeon duly qualified and trained in the polyactic interference screw use (knee ligament reconstruction) and who has read these instructions for use. The entire packaging of the device ensures a physical barrier essential for maintaining its sterile state. The polyactic interference screw is thus sterile unless the package has been opened, damaged or otherwise contaminated. Before opening the package, the surgeon and/or his team must verify that there is no sign of leakage, hole, tear or defective sealing of the package that would require disposal of the device. Any contamination of the surgical equipment should be minimized to prevent postoperative infection. Before implanting the device, the surgeon must ensure that active or recent infections have been treated properly.

Please refer to Surgical Technique documentation available from the manufacturer.

Any product that has been in contact with the patient, tissues or organic fluids must be considered as clinical waste and must be disposed of in accordance with local regulations for the disposal of clinical waste.

RECOMMENDATIONS FOR USE

preoperative recommendations

The surgeon is responsible for following preoperative and postoperative procedures and for knowing and mastering the surgical technique used. The medical device and the associated instrumentation must be selected prior to the surgical procedure. The patient must be informed of the preoperative and postoperative procedures as well as of the contraindications and adverse effects that may be observed.

Recommendations for use during the intervention

The polyactic interference screw must be removed from its packaging and handled with sterile gloves or instruments. Installation instructions appropriate for the selected surgical technique must be used for implant site preparation. The screwdriver Ref16NT001 designed specifically for compatibility with the screw internal imprint must be used to finalize the placement of the device and ensure fixation of the transplant. The screwdriver is not fully inserted into the screw or the screwdriver and screw are not properly aligned, both may be damaged. For femoral fixation, tighten the screw at the cortical bone.

For tibial fixation, insert the tibial screw at the anterior cortical bone in the distal part of the tibial tunnel.

Postoperative recommendations

The fixation provided by the polyactic interference screw is temporary until complete bone healing is achieved. Therefore, the device should not be subjected to loading or other stress for the duration of the time recommended by the surgeon during postoperative follow-up.

Use of the device to remove the polyactic interference screw should take into account the potential risk to the patient in the event of a second surgical procedure. Removal of the device should be followed by an appropriate postoperative procedure.

DEVICE RELATED INCIDENTS

Any serious incident in relation to the use of the device should be reported to the manufacturer and to the competent authority of the Member State where the event took place.

The manufacturer must be contacted at the following email address: materiovigilance@biomatlante.com

The contact points of the competent authorities are available on the website of the European Commission: <https://ec.europa.eu/tools/eudamed>.

ADDITIONAL INFORMATION

The basic UDI-DI number assigned to the polyactic interference screw is: 37050518213040133

The databases below can be consulted in order to obtain additional information concerning the device:

- Summary of safety and clinical performance: <https://www.admedsol.com/information/>
- Instructions for use: <https://www.admedsol.com/information/>

REF	Catalogue number	MD	Medical Device	MR	MR safe
LOT	Batch code	STERILE N	Sterilized using irradiation		Caution
	Date of manufacture FR: Country of manufacture (France)		Do not re-sterilize		CE mark
	Use-by date		Do not re-use		Consult instructions for use or consult electronic instructions for use
	Unique Device Identifier		Do not use if package is damaged and consult instructions for use		Patient information website
	Manufacturer		Upper limit of temperature 32°C (90°F)		Double sterile barrier system
	Distributor		Keep dry		

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COMPOSITION

Composant	Proportion
Acide poly(L-Lactide-co-D-Lactide)	100%

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