

# Biddu

## High Tibial Osteotomy System




In the treatment of lower extremity malalignment and osteoarthritis, high tibial osteotomy is a procedure that allows for the realignment of the lower extremity and thereby either shifts the weight-bearing access to the healthier portion of the knee or corrects the alignment deformity.

This procedure is performed by cutting the tibia and with the use of an opening jack allows the surgeon to dial in the exact correction desired to correct the angular deformity or shift the weight bearing axis to a desired position in the knee.

An anatomically precontoured low-profile A/P wedge plate and low-profile head screws is applied to fixate the osteotomy.

### Instruments Set



### Product List

BAK-7094.075S	BIDDU Tibial A / P Sloped Wedge Osteotomy Plate, 7.5mm, SS, Sterile
BAK-7094.100S	BIDDU Tibial A / P Sloped Wedge Osteotomy Plate, 10.0mm, SS, Sterile
BAK-7094.125S	BIDDU Tibial A / P Sloped Wedge Osteotomy Plate, 12.5mm, SS, Sterile
BAK-7094.150S	BIDDU Tibial A / P Sloped Wedge Osteotomy Plate, 15.0mm, SS, Sterile
BAK-7181.45.28	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 28mm, SS, Sterile
BAK-7181.45.30	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 30mm, SS, Sterile
BAK-7181.45.32	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 32mm, SS, Sterile
BAK-7181.45.34	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 34mm, SS, Sterile
BAK-7181.45.36	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 36mm, SS, Sterile
BAK-7181.45.38	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 38mm, SS, Sterile
BAK-7181.45.40	BIDDU LPH Cortex Screw, Self tapping, Dia. 4.5mm, Length 40mm, SS, Sterile
BAK-7182.65.45	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 45mm, SS, Sterile
BAK-7182.65.50	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 50mm, SS, Sterile
BAK-7182.65.55	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 55mm, SS, Sterile
BAK-7182.65.60	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 60mm, SS, Sterile
BAK-7182.65.65	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 65mm, SS, Sterile
BAK-7182.65.70	BIDDU LPH Cancellous Screw, Fully Threaded, Dia. 6.5mm, Length 70mm, SS, Sterile