

SLAP & Bankart Repair

MINI-VIM[®] Ligament Anchor (Titanium)

The MINI-VIM[®] Ligament Anchors are titanium anchors available in 2.8mm and 3.5mm diameters and come preloaded with USP #2 BioFiber[®] suture. In most cases, the implant can be manually inserted using the handled version. A drill is available for use in hard bone if desired. For arthroscopic applications, the MINI-VIM[®] anchor can be inserted through the small diameter MINI-VIM[®] spear eliminating the need for a cannula. This is ideal for SLAP and subscapularis repairs where a small stab incision and percutaneous delivery of the implant is preferred. MINI-VIM[®] Ligament Anchor is screw-in anchor.

Key Features:

- Tremendous pull-out strength.
- Self-drilling. A drill is available for use in hard bone if desired.
- Preloaded with One USP #2 BioFiber[®] in a single-use sterile pack.
- Also available in needed version.

MINI-VIM[®] Ligament Anchor

BAS-9001.28F	MINI-VIM [®] Ligament Anchor, Titanium, loaded with 1 pc Fiber; Dia. 2.8mm
BAS-9040.28F	MINI-VIM [®] Ligament Anchor, Titanium, loaded with 1 pc Fiber with needle, Dia. 2.8mm
BAS-9031.35F	MINI-VIM [®] Ligament Anchor, Titanium, loaded with 1 pc Fiber, Dia. 3.5mm
BAS-9083.35F	MINI-VIM [®] Ligament Anchor, Titanium, loaded with 1 pc Fiber with needle, Dia. 3.5mm

Related Instrumentation

BAS-9012	MINI-VIM [®] Spear (Two-pronged) for BAS-9001.28F & BAS-9031.35F
BAS-9109	MINI-VIM [®] Spear (Fishmouth) for BAS-9001.28F & BAS-9031.35F
BAS-9013	Trocar for BAS-9012 & BAS-9109
BAS-9026	Drill bit, dia. 2.2mm for MINI-VIM [®] (BAS-9001.28F)
BAS-9114	Drill bit, dia. 2.8mm for MINI-VIM [®] (BAS-9031.35F)



SLAP & Bankart Repair

MINI-VIM PK[®] Ligament Anchor (PEEK)

The MINI-VIM PK[®] Ligament Anchors are available in a variety of sizing and suture options. The MINI-VIM PK[®] is a nonabsorbable PEEK suture anchor with a material eyelet which provides superior abrasion resistance due to PEEK's low coefficient of friction. Anchor inserting the anchor through the same spear into the drilled hole. MINI-VIM PK[®] Ligament Anchor is Tap-in/Bang-on anchor



MINI-VIM PK[®] Ligament Anchor

BAS-9086.28F	MINI-VIM PK [®] Ligament Anchor, PEEK, loaded with 1 pc Fiber, Dia 2.8mm
BAS-9087.28F	MINI-VIM PK [®] Ligament Anchor, PEEK, Knotless, Dia. 2.8mm
BAS-9088.30F	MINI-VIM II PK [®] Ligament Anchor, PEEK, loaded with 2 pc Fiber, Dia. 3.0mm
BAS-9089.35F	MINI-VIM II PK [®] Ligament Anchor, PEEK, loaded with 2 pc Fiber, Dia. 3.5mm
BAS-9192.35F	MINI-VIM PK [®] Ligament Anchor, PEEK, Knotless, Dia. 3.5mm

Related Instrumentation

BAS-9012	Spear (Two-pronged) for MINI-VIM [®] BAS-9001.28F
BAS-9109	Spear (Fishmouth) for MINI-VIM [®] BAS-9001.28F
BAS-9013	Trocar for BAS-9012 & BAS-9109
BAS-9110	Spear (Two-pronged) for MINI-VIM PK [®] BAS-9088.30F & BAS-9089.35F
BAS-9097	Spear (Fishmouth) for MINI-VIM PK [®] BAS-9088.30F & BAS-9089.35F
BAS-9098	Trocar for BAS-9097 & BAS-9110
BAS-9100	Drill bit, Dia. 2.4mm for MINI-VIM PK [®] (BAS-9086.28F)
BAS-9111	Drill bit, Dia. 2.5mm for MINI-VIM PK [®] , Knotless (BAS-9087.28F)
BAS-9101	Drill bit, Dia. 2.6mm for MINI-VIM II PK [®] (BAS-9088.30F)
BAS-9102	Drill bit, Dia. 3.0mm for MINI-VIM II PK [®] (BAS-9089.35F)
BAS-9127	Sleeve for MINI-VIM PK [®] , Knotless Anchor, Dia. 2.8mm