

C-Thru

Disposable Cannula-Trocar Set



- Threaded translucent twist-in design
- Rounded distal tip
- Allows direct arthroscopic visualization of instruments and suture passing through the cannula
- Flexible option easily conforms to large or curved instruments
- Supplied with a disposable obturator
- Ratcheting stopcock for controlled outflow
- Latex-free design
- Available in pre-sterilized packaging
- Reusable Obturators available in stainless steel for use with disposable cannulas

Cat. no. Product

BAS-9071	C-Thru Disposable cannula-trocar set, Dia. 5.5mm, Len. 90mm, Sterile
BAS-9072	C-Thru Disposable cannula-trocar set, Dia. 6.5mm, Len. 70mm, Sterile
BAS-9073	C-Thru Disposable cannula-trocar set, Dia. 8.5mm, Len. 70mm, Sterile
BAS-9074	C-Thru Disposable cannula-trocar set, Dia. 8.5mm, Len. 90mm, Sterile
BAS-9004.08	Reusable Cannula Introducer, 8.0mm
BAS-9082	Reusable Cannula Introducer 5.5mm