

Flexi Cannulated Reamers

For Knee Arthroscopy





Provides flexibility for anatomic cruciate ligament reconstructions



The unique design of the Flexi Cannulated Reamers provides an accurate and reproducible system for anatomic cruciate ligament reconstruction



The puzzle cut design offers excellent flexibility

Key Features:

- Reaches more inferior and anterior sites on the lateral condyle than conventional reamers
- Low Profile head enables precise intra-articular placement and safe tunnel
- No need to hyperflex the knee while drilling through medial portal
- · Allows anatomic placement of femoral tunnel
- Reamers are available in 5.0mm to 12.0mm diameters with 0.5mm increments
- Cannulated design to be used with 2.4mm flexible guide pin
- 4.7mm reamer is available for use with ONButton CL® and Buttonfix® fixation techniques







Ordering Information

Flexi Cannulated Reamer Set [12 pcs], Dia. 4.7mm to 10.0mm

*Flexi Cannulated Reamers are available in set only and available for Indian market only .







BT/AR/20, Rev. 01