

SOFTFIX-PK[®] Interference screw (PEEK)



- SOFTFIX-PK[®] Interference Screws provide strong mechanical fixation for both bone patellar tendon-bone (BTB) and soft tissue grafts in ACL and PCL reconstruction.
- SOFTFIX-PK[®] Interference Screws feature a thread pattern that allows for a simple surgical technique, with minimal tunnel preparation.
- Non-absorbable SOFTFIX-PK[®] interference screw incorporates PEEK-VESTAKEEP[®] polymer from Evonik[®].
- PEEK-VESTAKEEP[®] material exhibits an ideal combination of strength, stiffness and toughness which surpasses absorbable materials.
- PEEK (Poly-Ether-Ether-Ketone) is a thermoplastic material with excellent biocompatibility and bio stability characteristics as per ASTM F2026.
- SOFTFIX-PK[®] interference screws offer the revision ability of an absorbable screw.
- PEEK has inherent radiolucent characteristics that have widely gained confidence and acceptance in the medical community.

SOFTFIX-PK[®] Interference screw, PEEK

BAK-7164.15	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 5.0mm, Len. 15mm
BAK-7164.20	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 5.0mm, Len. 20mm
BAK-7165.23	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 6.0mm, Len. 23mm
BAK-7166.23	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 7.0mm, Len. 23mm
BAK-7166.28	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 7.0mm, Len. 28mm
BAK-7167.23	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 8.0mm, Len. 23mm
BAK-7167.28	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 8.0mm, Len. 28mm
BAK-7168.23	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 9.0mm, Len. 23mm
BAK-7168.28	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 9.0mm, Len. 28mm
BAK-7168.35	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 9.0mm, Len. 35mm
BAK-7169.28	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 10.0mm, Len. 28mm
BAK-7169.35	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 10.0mm, Len. 35mm
BAK-7170.35	SOFTFIX-PK [®] Interference screw, cannulated, roundhead, Dia. 11.0mm, Len. 35mm

Related Instrumentation

BAK-7191	Guide Wire for SOFTFIX-PK [®] Interference Screws Dia. 1.2mm; Len. 300mm (SS)
BAK-7180	Guide Wire for SOFTFIX-PK [®] Interference Screws Dia. 1.2mm; Len. 300mm (Nitinol)
BAK-7161	Screwdriver for SOFTFIX-PK [®] PEEK Interference Screws Cannulated

SOFTFIX[®] Interference screw (Titanium)



- Specifically designed for endoscopic femoral fixation where graft protection is needed.
- Polished, fully radiused threads provide maximum graft protection.
- Wide-bore cannulation for easy guide wire removal.
- The Taper-Tip style offers sharper initiating threads and a more pronounced taper for faster, easier insertion.
- Designed to reduce the possibility of bone plug migration during placement.
- Round head reduces possibility of graft damage during or after implantation.
- Universal screw for either tibial or femoral applications. Standard 4.5mm hex for easy insertion and removal.
- Cannulated for use with a guide wire.
- Material - Titanium alloy.

SOFTFIX[®] Interference screw, Titanium

BAK-7128.20S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 7.0mm x Len. 20mm
BAK-7128.25S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 7.0mm x Len. 25mm
BAK-7128.30S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 7.0mm x Len. 30mm
BAK-7101.20S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 8.0mm x Len. 20mm
BAK-7101.25S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 8.0mm x Len. 25mm
BAK-7101.30S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 8.0mm x Len. 30mm
BAK-7102.25S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 9.0mm x Len. 25mm
BAK-7102.30S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 9.0mm x Len. 30mm
BAK-7114.25S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 10.0mm x Len. 25mm
BAK-7114.30S	SOFTFIX [®] Interference screw, cannulated, roundhead, Dia. 10.0mm x Len. 30mm

Related Instrumentation

BAK-7120	Guide Wire for SOFTFIX [®] Interference Screws Dia. 2.0mm; Len. 300mm (SS)
BAK-7183	Guide Wire for SOFTFIX [®] Interference Screws Dia. 2.0mm; Len. 300mm (Nitinol)
BAK-7063	Screwdriver for SOFTFIX [®] Titanium Interference Screws Cannulated

