

## **Buttonfix**<sup>®</sup>

### Fixation Button with Adjustable Loop

One of the strongest adjustable suspensory fixation devices for cruciate reconstruction



# Buttonfix<sup>®</sup>

## Fixation Button with Adjustable Loop

### Adjustable Suture

- Enables calibrating the loop to its optimum size.
- Allows the surgeon to maximize the amount of tendon in the femoral tunnel, thereby optimizing the graft-tunnel healing interface and the chance of a successful outcome.

### Green suture for flipping

- To flip the button over the lateral cortex of femur.

### Titanium Button

- Strong with low profile.
- Compatible with common extension devices.

### UHMWPE Suture Loop

- Provide high strength to the loop.
- As strong as Closed loop.

### Graft Support Frame

- A larger and thicker surface contact with the graft.
- Protects the graft from damage and any abrasion.

### White suture for pulling lead

- To pull the graft from tibial and femoral tunnel.

### Ordering Information

BAK-7195 BUTTONFIX<sup>®</sup> Fixation Button with Adjustable Loop



**BIOTEK - Chetan Meditech Pvt. Ltd.**  
Plot No: MD 04, Near Teva Company, Sanand GIDC Gate No. 2,  
Charal Gaam, Ahmedabad-382110. Gujarat, INDIA.  
E: [contact@biotekortho.com](mailto:contact@biotekortho.com); W: [www.biotekortho.com](http://www.biotekortho.com)

• ISO 13485:2016  
• CE Certified

